

RADIOPAQUE GLASS IONOMER LUTING CEMENT (HUMID STABLE)

ISO 9917-1:2007 Glass Polyalkenoate Class 4.2a & 4.2b

DIRECTIONS FOR USE

Adheres chemically to tooth substance and non-precious alloys, which gives maximum secure retention. Excellent biocompatibility. Contains fluoride ions. Radiopaque.

INDICATIONS:

Use as liner or base.

For cementing crowns, bridges, inlays, and onlays.

For cementing orthodontic brackets and bands.

MIXING:

The recommended powder: liquid ratio is 2:1 (mm) at 22-24°C and 40-60% R.H. The use of a clean dry, polished glass slab or paper pad and a stainless steel "Silicate" spatula is recommended. Invert bottle to "fluff" powder for accurate dispensing. Measure ONE level scoop of powder onto a glass slab or paper pad. Dispense ONE "bubble-free" drop of liquid onto the glass slab or paper pad. Incorporate half the powder into the liquid and mix for 15 seconds, then add the remaining powder and mix to a uniform syrup-like consistency. Total mixing time is 30 seconds.

DO NOT ADD POWDER IN SMALL INCREMENTS.

WORKING TIME: 2 minutes 30 seconds (minimum) from start of mix at 22-24°C and 40-60% R. H.

CLINICAL SETTING TIME: A maximum of 6.5 minutes from start of mix.

APPLICATION:

Apply the mixed cement immediately on completion of mix. Excess cement may be removed immediately or after cement has set. DO NOT attempt to remove excess cement whilst cement is in the "Rubbery" phase of setting (between 2.5-6 minutes from start of mix) as adhesion may be impaired. Wipe the filling instrument/spatula clean immediately.

STORAGE:

Store powder in a cool, dry place. The powder is protected against moisture uptake by a desiccant capsule, since moisture ingress can adversely affect the shelf life. Always replace cap immediately after use. If desiccant becomes pink, discard the product.

WARNING:

Avoid contact of liquid and powder with eyes. In the event of contact, rinse with copious quantities of water and seek medical assistance.

BATCH CODE:



The batch code gives an open date of manufacture in month, year, day format with a numerical suffix to uniquely identify the batch of material. Please quote this batch number in all correspondence.



The expiry date is shown in year, month format. Do not use the product after this date.

TESTING:

For test purposes, the powder: liquid ratio is: 1.00g: 0.50g at 22-24°C. and 40-60% R.H.

ISO 9917 Net Setting Time: A maximum of 5 minutes 30 seconds from end of mix at 36-38°C.

SPECIALLY FORMULATED FOR USE IN DENTISTRY

REF:	259XPL15	15g Powder / 8ml Liquid
	259X15	15g Powder / 10ml Liquid
	259XPL	30g Powder / 15ml Liquid



Advanced Healthcare Ltd., Tonbridge, Kent, TN11 8JU, UK

Tel: +44 1892 870500

Fax: +44 1892 870482

Advanced Healthcare Ltd operate a policy of continuous monitoring and improvement of our products. If you have any comments about this product, please contact us at the above address stating the batch number shown on the packaging.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a dental professional.