

REINFORCED ZINC OXIDE EUGENOL FILLING MATERIAL

ISO 3107 Type II

DIRECTIONS FOR USE

The Zinc Oxide Eugenol filling material is reinforced with an acrylic resin to give excellent strength in order to withstand condensation forces, as a lining material and mastication forces as a temporary restorative.

INDICATIONS:

Base or lining under restorative materials (amalgam, silicate, silicophosphate, glass ionomer).

Temporary restorative.

CONTRA-INDICATIONS:

As with other Eugenol containing materials, the polymerisation of acrylic direct filling materials (acrylic resins and composite resins) and temporary crown and bridge materials is inhibited and use with these materials must be avoided.

MIXING:

Recommended ambient conditions for mixing 22-24°C, 48-52% RH. Mix ratio 4:1 m:m.

Dispense powder and liquid onto a glass slab. Divide the powder into five parts. Using a steel spatula, mix in increments at 15-second intervals. Spatulate completed mix for 15 seconds. Total mixing time 1½ minutes.

WORKING TIME: 2 minutes from end of mix.

SETTING TIME: A maximum of 6 minutes at 36-38°C from start of mix.

WARNING:

Liquid contains Eugenol. IRRITANT. Irritating to skin. Avoid contact with skin and eyes.

STORAGE:

Store out of sunlight in a cool dry place 5-25°C.

BATCH CODE:



The batch code gives an open date of manufacture in month, year, day format with a numerical suffix to uniquely identify the batch of material. Please quote this batch number in all correspondence.



The expiry date is shown in year, month format. Do not use the product after this date.

SPECIALLY FORMULATED FOR USE IN DENTISTRY

| | | |
|------|--------|--------------------------|
| REF: | 046XPL | 40g Powder / 15ml Liquid |
| | 046X1 | 40g Powder |
| | 127X3 | 15ml Liquid |



Advanced Healthcare Ltd., Tonbridge, Kent, TN11 8JU, UK

Tel: +44 1892 870500

Fax: +44 1892 870482

Advanced Healthcare Ltd operate a policy of continuous monitoring and improvement of our products. If you have any comments about this product, please contact us at the above address stating the batch number shown on the packaging.

Caution: U.S. Federal Law restricts this device to sale by or on the order of a dental professional.

2017-02

AP8164