

AMALGOMER™ CR

Technology

IMPROVED

Glass Ionomer Ceramic Reinforced Posterior Restorative (Powder Liquid Type)

ISO 9917-1 GLASS POLYALKENOATE

DIRECTIONS FOR USE

FEATURES:

Improved Amalgomer CR is an advanced posterior restorative material that contains fluoride, is radio opaque and has exceptional strength. It bonds chemically to enamel, dentine and cementum. Passes the tests in the Amalgam normative standard ISO 24234:2004.

Now translucent.

INDICATIONS:

- Class I and II cavities in deciduous teeth.
- Repair of amalgam restored teeth when either tooth or restoration has fractured.
- Class I and II cavities in selected permanent teeth.
- Base under amalgam and posterior composite restorations.
- In all classes of cavity where radiopacity is a prime requirement.
- As a core build up under crown.
- On the root surfaces for locating overdentures.
- Long term temporary replacement for cusp(s).
- Repair to crown margins.

CAVITY PREPARATION:

Use minimal tooth reduction wherever possible. Calcium hydroxide liners need only be used in deep cavities. Tooth Cleanser (25% polyacrylic acid solution) is applied to areas of unreduced tooth, to which adhesion is required, using a pledget of cotton wool for a maximum of 30 seconds. Wash with water and dry with oil-free air. If freshly cut enamel or dentine is contaminated with saliva, Tooth Cleanser should be applied but only for 10 seconds before washing is started.

MIXING:

Powder liquid ratio 3.6:1 m:m at 22-24°C and 40-60%RH).

Use a glass block for best results and a stainless steel "Silicate" spatula. Incorporate half the powder into the liquid as quickly as possible (5-10 seconds) and then add the remainder and spatulate to a THICK putty-like consistency. Total mixing time 30 seconds. DO NOT ADD POWDER IN SMALL INCREMENTS.

PLACEMENT:

Apply to cavity using normal instruments (non-stick aluminium instruments are particularly suitable). If problems of cement adhering to the steel instrument occur, dip the clean instrument either into the powder or methylated spirit. If stainless steel matrices are used, they should be lightly coated with petroleum jelly.

Working time: A minimum of 1 minute 30 seconds from start of mix at 22-24°C.

Clinical Set Time: A maximum of 3 minutes 30 seconds from end of mix.

FINISHING:

Best results are obtained by polishing with abrasive discs and stones at a subsequent visit using water spray lubrication. At 7 minutes after placement, the material is hard enough to finish using abrasive discs and stones, but petroleum jelly must be used as a lubricant to prevent excess heat and desiccation of the cement. After finishing, the surface should be coated with varnish for protection from moisture.

WARNING:

Avoid contact of cleanser liquid and powder with the eyes. In case of such contact, wash thoroughly with large quantities of water and obtain medical assistance.

STORAGE:

The powder is protected against moisture uptake by a desiccant capsule, since moisture ingress can adversely affect the shelf life. Always replace caps immediately after use. Store in a cool dry place (5-25°C).

BATCH CODE:



The batch code gives an open date of manufacture in month, year, day format with a numerical suffix to uniquely identify the batch of material. Please quote this batch number in all correspondence.



The expiry date is shown in year, month format. Do not use the product after this date.

† For most purposes the powder/liquid ratio is 0.72g to 0.20g at 23±1°C and 50± 10% RH.

Also passes the tests in the ISO standard for Amalgam (ISO 24234:2004).

SPECIALLY FORMULATED FOR USE IN DENTISTRY

REF: 12g Powder / 5ml Liquid

Amalgomer CR also passes the tests in the ISO standard for Amalgam (ISO 24234:2004)

Specially formulated for use in dentistry.



Advanced Healthcare Ltd., Tonbridge, Kent, TN11 8JU, UK
Tel: +44 1892 870500 Fax: +44 1892 870482
E-mail: sales@ahl.uk.com Web: www.ahl.uk.com



Advanced Healthcare Ltd. operates a policy of continuous monitoring and improvement of our products. If you have any comments about this product, please contact us at the above address stating the batch number shown on the packaging.

2017-06
AP7021

THE STRENGTH
BEHIND YOUR SMILE