



FEEDBACK
www.ahl.uk.com/feedback

AHL

reliaFlow Bulk

Bulk Flow Composite Restorative

Complies With ISO 4049

DIRECTIONS FOR USE

FEATURES:

An advanced flowable composite predominantly designed as a base under class I & II restorations. Can be placed in 4mm layers. Self-levelling action allows easy adaptation. High filler content to reduce polymerisation shrinkage. Radiopaque. Contains fluoride.

INTENDED PURPOSE:

Base or liner under other restorative materials

INTENDED PATIENT POPULATION:

From child to geriatrics

INTENDED USER:

This product has been formulated for use in dentistry and is intended for use by dental professionals only.

CLINICAL BENEFIT:

To restore the function of the teeth and help maintain the integrity of the remaining tooth structure.

INDICATIONS FOR USE:

- Base in class I & II direct restorations.
- Liner under direct restorative materials.
- Restorations of small posterior cavities (including occlusal restorations).

CONTRA-INDICATIONS:

Do not use this product for any purposes other than specifically outlined under INDICATIONS in these instructions for use

CONTENTS OF PACK:

Syringe 2g, TIPS x 5, Instructions for use

PRECAUTIONS AND WARNINGS:

- Do not expose patients or users known to be allergic to this type of material and/or methacrylate monomers.
- If any eruption or skin inflammation occurs on a patient when using this product, immediately discontinue use and have the patient seek medical advice.
- Avoid contacting oral mucosa, skin or eyes. In case of contact, rinse immediately with copious amounts of water and seek medical advice.
- Monomers may pass through dental gloves; avoid material contact for maximum protection.
- Wear UV protective glasses when light curing.
- When storing these materials in a refrigerator, allow them to reach room temperature before use

PROCEDURE

(1) CAVITY PREPARATION:

- Thoroughly clean the tooth surface to remove plaque.
- Prepare the cavity in the conventional way.
- If the pulp is exposed or is close to the restoration, then a calcium hydroxide liner should be used. (N.B. Do not use a eugenol containing liner).
- Apply the bonding agent according to the manufacturer's directions.
- AHL Generic Etch Gel & Bonding Agent is recommended (code: 071X1 & 250X1) for use in cavities.

(2) PLACEMENT:

- Fit dispensing tip and dispense the required amount of material from the syringe directly in to the bottom of the cavity.
- Take care not to include air bubbles when dispensing material.
- Always replace the syringe cap after use.
- Place in a layer of up to 4mm depth and cure with a dental light cure unit.

(3) CURING

- Select the light curing time in the table below.

Shade	Depth of cure	LED curing unit	Halogen lamp curing unit
		Wavelength 440-490 nm	Wavelength 400-500 nm
		Light Intensity $\geq 1000 \text{ mW/cm}^2$	Light Intensity $\geq 500 \text{ mW/cm}^2$
Universal	4 mm	10 seconds	20 seconds
Dentine		20 seconds	40 seconds

- Use a dental light-curing unit having specifications equivalent to those described above.
- When using a unit having performance other than these recommendations, follow the instructions for that unit

(4) FINISHING:

- After curing, cap with AHL universal composite (UL316)

SINGLE USE

Discard the tip after use.
This must not be re-used on another patient for hygienic reasons.

STORAGE:



Store in a cool, dry place (1-30°C).
Always replace cap immediately after use.
Keep away from direct sunlight

EXPIRY:



The expiry date is shown in year, month format.
Do not use the product after this date.

DISPOSAL:

Dispose of the contents and containers in accordance with relevant local and national requirements.

POSSIBLE SIDE EFFECTS / RESIDUAL RISKS:

- This product contains substances that may cause an allergic reaction.
- Restorations have the potential to fracture depending on patient habits.
- Restorations have the potential to fall out depending on patient habits.

BATCH CODE:

LOT

The batch code gives an open date of manufacture in month, year, day format with a numerical suffix to uniquely identify the batch of material. Please quote this batch number in all correspondence.

DEVICE CODES:

REF

AH0920

-

2g Syringe Universal

AH0921

-

2g Syringe Dentine

COMPOSITION:

Composition	% by weight
Filler	55 – 85
Methacrylate resins	10 – 45
Initiators	<1
Stabilisers	<1
Pigments	<1
Other	<1



AHL operate a policy of continuing surveillance & monitoring of our products. If you experience any incidents relating to the use of this product, please immediately contact us at the below address stating the batch number shown on the packaging. If you experience any serious incident relating to the use of this product, please immediately contact AHL at the below address and the competent authority of the territory you are in.

A summary of safety & clinical performance (SSCP) is available via the EUDAMED database.
<https://ec.europa.eu/tools/eudamed>

Caution: U.S. Federal Law restricts this device to sale by or on the order of a dental professional.



Advanced Healthcare Ltd., Tonbridge, Kent, TN11 8JU, UK
Tel: +44 1892 870500 Email: sales@ahl.uk.com

EU REP

Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatara, BKR 4013 Malta.

2025-05

AP9386/3