



# AHL Generic



## TEMPORARY LUTING CEMENT

### DIRECTIONS FOR USE

#### FEATURES:

A ready to use temporary cement which fixes dental crowns and bridges. Cement sets in contact with water or saliva. It does not irritate the gums, it is easy to work with and can be easily removed without damage to the bridge or crown. It is not a permanent filling material and recommended use is up to a maximum of 4 weeks only.

#### INTENDED PURPOSE:

Temporarily attach dental prostheses.

#### INTENDED PATIENT POPULATION:

From child to geriatrics

#### INTENDED USER:

This product has been formulated for use in dentistry and can be use by dental professionals or the general public.

#### CLINICAL BENEFIT:

To temporarily restore the function of the teeth and help maintain the integrity of the remaining tooth structure.

#### INDICATIONS FOR USE:

- A temporary cement for fixing dental bridges or crowns.

#### CONTRA-INDICATIONS:

- Do not use for more than 4 weeks

#### CONTENTS OF PACK:

262X28G 4 x 7g tubes, instructions for use  
262X8.4 3 x 2.8g syringes, 9 x tips, instructions for use  
224X2 1 x 7g tube, instructions for use

#### PRECAUTIONS AND WARNINGS:

- Do not expose users known to be allergic to this type of material.
- To avoid infection and/or serious injury, this product should not be used if there is throbbing pain, ulceration or swelling in the affected area.
- If any of these symptoms exist, consult your dentist immediately.
- This product should not be placed on open wounds.
- READ ALL OF THE DIRECTIONS CAREFULLY BEFORE USE!
- Wash and dry your hands.
- Do not get unused cement wet.
- Use only what is needed to repair the dental emergency.
- KEEP OUT OF REACH OF CHILDREN.
- DO NOT use product for any purpose other than indicated.
- You should consult a dentist within 24 hours of applying this product.
- Temporary filling material is only used as a temporary solution.

#### PROCEDURE

##### (1) PREPARATION:

- Clean out the crown or bridge (if necessary) and leave moist but not wet.
- Clean the area in the mouth (tooth stump) where the crown is to be replaced.
- Rinse the mouth with water (do not dry the stump).

##### (2) PLACEMENT:

- Apply a small amount of the cement around the inner edge of the crown or cap.
- Do not fill the crown.
- For bridges, apply a small amount of the cement around the inner edge of each crown of the bridge.
- Fit back in the mouth and firmly press down. Then carefully bite into place. You will know when the crown or bridge is correctly re-seated as the teeth will close back together again comfortably.
- Press the crown or bridge in place for 3 minutes. (If you are not satisfied with the fit, you have up to 5 minutes to remove the bridge or crown and repeat the whole process).
- Remove any excess cement from around the base of the crown.
- Do not chew for the next 2 hours.

#### STORAGE:



Store in a cool, dry place (5-25°C).  
Always replace cap immediately after use.

#### EXPIRY:



The expiry date is shown in year, month format. Do not use the product after this date.

#### DISPOSAL:

Dispose of the contents and containers in accordance with relevant local and national requirements.

#### POSSIBLE SIDE EFFECTS / RESIDUAL RISKS:

- This product contains substances that may cause an allergic reaction.

#### BATCH CODE:



The batch code gives an open date of manufacture in month, year, day format with a numerical suffix to uniquely identify the batch of material. Please quote this batch number in all correspondence.

#### DEVICE CODES:



262X28G 4 x 7g Tube  
262X8.4 3 x 2.8g Syringe Kit  
224X2 1 x 7g Tube

#### COMPOSITION:

Composition	% by weight
Calcium Sulphate	20 – 70
Barium Sulphate	0 – 10
Glycol Ether	10 – 30
Zinc Oxide	5 – 25
Acrylic Polymer	0 – 5
Potassium Sulphate	0 – 20

AHL operate a policy of continuing surveillance & monitoring of our products. If you experience any incidents relating to the use of this product, please immediately contact us at the below address stating the batch number shown on the packaging. If you experience any serious incident relating to the use of this product, please immediately contact AHL at the below address and the competent authority of the territory you are in.

A summary of safety & clinical performance (SSCP) is available via the EUDAMED database. <https://ec.europa.eu/tools/eudamed>

Caution: U.S. Federal Law restricts this device to sale by or on the order of a dental professional.



Advanced Healthcare Ltd., Tonbridge, Kent, TN11 8JU, UK  
Tel: +44 1892 870500 Email: [sales@ahl.uk.com](mailto:sales@ahl.uk.com)



Advena Ltd. Tower Business Centre, 2nd Flr., Tower Street, Swatar, BKR 4013 Malta.

2025-05

AP8060/4

