

AHL

tubliCEM

Zinc Oxide – Eugenol Root Canal Sealer

Product Description

A paste-paste zinc oxide–eugenol cement for use with obturation points. The 2 pastes mix easily to a smooth creamy consistency with long working time and setting time. The cement has moderate flow characteristics, is radiopaque and is low in film thickness (about 6µm).

Features

- Excellent flow characteristics
- Long working time during preparation
- Zinc Oxide Eugenol formula
- Radiopaque
- Tube – Tube presentation

Benefits

- Easy to apply, flows well into canal
- Aids material manipulation outside mouth
- Low leakage and low dimensional change
- Easily distinguished from surrounding tooth structure
- Easy to prepare mix and convenient to use.

Indications

- Permanent obturation of the root canal space with the aid of obturation points.

Product Presentation: 8g Base Paste / 3.5g Catalyst Paste (in tubes)

Shelf life: 3 years from the date of manufacture when supplied in AHL's standard packaging.