





- Adhesive to enamel & dentine.
- ✓ Proven formula
- Wide range of uses
- ✓ Radiopaque

ZCem Phos Zinc Phosphate Cement

ZCem Phos combines tried and tested technology with modern processing techniques to produce a dependable and durable cement.

Although first introduced in the 1890's zinc phosphate cements still play an important role in modern dentistry.

Indications

Cementation of crowns and bridges Cementation of inlays and onlays Cementation of Orthodontic bands

Features

Adhesive to enamel and dentine Proven formula Radiopaque

Advantages

Wide range of uses Shows on X-rays

Available in	Code
Powder/Liquid 90gm/30ml	AH0600

Property	ISO 9917-1:2007 Class 4.2.a	Typical Properties
Working Time	-	3 mins
Net Setting Time	2 min 30 secs – 8 mins	5 mins 10 secs
Compressive Strength	50MPa minimum	148MPa
Film Thickness	25µm Maximum	21μm

Your global partner

Delivering quality for more than 30 years

Advanced Healthcare Ltd.

Units 2-4 Leavers Estate, Chiddingstone Causeway, Tonbridge, Kent TN11 8JU UK Tel: +44(0)1892 870500 Fax: +44(0)1892 870482 Email: sales@ahl.uk.com

www.ahl.uk.com

© 2023 Advanced Healthcare Ltd. E&OE.

Authorised distributor