



- ✓ Mixes easily
- ✓ Good handling
- ✓ Radiopaque
- ✓ Temporary restorative

ZOECem Reinforced Zinc Oxide Eugenol Cement

ZOECem features a reinforcing polymer incorporated into the powder, giving the cement the strength to resist condensation forces and to ensure adequate life when used as a temporary filling.

By incorporating the reinforcing agent in the powder instead of the liquid, the mixing properties are not affected. Although the cement is reinforced, it can be easily removed when used in a temporary situation.

Indications

Temporary filling material
Base or lining under restorative materials (amalgam, silicate, silicophosphate, glass ionomer)

Features

Reinforcing component present in powder
Radiopaque
Temporary restorative.

Advantages

Good mixing properties
Shows on X-rays

Available in	Code
40gm Powder/15ml Liquid	AH0750

Property	ISO 3107:2011 Type II	Typical Properties
Working Time	-	4 mins
Net Setting Time	1 min 30 secs – 10 mins	4 mins 45 secs
Compressive Strength	5MPa minimum	38MPa
Disintegration	-	0.50%

Your global partner

Delivering quality for more than 30 years

Advanced Healthcare Ltd.

Units 2-4 Leavers Estate, Chiddingstone Causeway, Tonbridge, Kent TN11 8JU UK
Tel: +44(0)1892 870500 Fax: +44(0)1892 870482 Email: sales@ahl.uk.com

www.ahl.uk.com

Authorised distributor